

## Laboratory Refrigerator

- Exterior made of sheet metal with a white epoxy coating that gives it high resistance to impacts and daily use.
- Inner body in ABS food material.
- They have interior fluorescent lighting that increases the appeal and exposure of the products.
- Curved joints on the inside to facilitate cleaning.
- Equipped with a key lock
- Ecological R-134a refrigerant.
- Digital thermometer, temperature control and automatic defrost.
- Working voltage: 230 V / 50 Hz.



## Technical Specifications

Model	SSLR-20
Internal Volume	350 litres
Temperature Range	-18°C
Form Factor	Stand Alone
Shelves	07 Nos. Steel coated
Included Components	1 Freezer
Configuration	Single Door
Wheels / Legs	Adjustable Feet
Color	White
Digital Indicator and controller	YES
Colour	Grey/White
Item Dimensions LxWxH	24 x 27 x 70 Inches
Power Supply	230V ±10%, 50Hz

**NOTE: Technical specification & design are subject to change without prior notice due to continuous upgradation**

### SEMCO SCIENTIFIC PRIVATE LIMITED

9C/91, Vrindavan Yojna, Raebareli Road, Near NKM School, Lucknow – 226029, U.P.

info@semcoscientific.in / semcoscientific@gmail.com / www.semcoscientific.in

+91-9918227744 / +91-6388224590 / 0522-3569875