

Digital Orbital Shaker

- Ergonomically designed with minimal footprint
- Maintenance -free permanent (self) lubricated shield ball bearings for noiseless operation
- FRP exterior with anti-vibration mount for smooth running
- Precise shaking speed control with gentle start and stop cycle
- Optional timer function for speed



Technical Specifications

Model	SSDOS-26
Platform	26*26 inch
Flask capacity	250ml * 48 flasks
RPM	200 RPM
RPM Controller	Speed Regulator
Construction	Heavy gauge mild steel powder coated
Clamps	Spring-loaded stainless steel
Power supply	220 Volts 50Hz
Optional	▪ Digital timer ▪ Digital RPM display

SEMCO SCIENTIFIC PRIVATE LIMITED

Address: 9C/91, Vrindavan Yojna, Raebareli Road, Near NKM School, Lucknow – 226029, U.P.
Email Id: info@semcoscientific.in / semcoscientific@gmail.com / **Website:** www.semcoscientific.in
Contact No.: +91-9918227744 / +91-6388224590 / 0522-3569875